Dentoclic

Glass fiber

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EN
Glass fiber posts for the reconstruction of endodontically treated teeth.

COMPOSITION (% EN PESO):

Gutta-percha 80%. Resina epoxi 20%.

INDICATIONS :

Prosthetic reconstruction of a tooth following root treatment.

PROPERTIES :

Esthetics. Reduces the risk of root fracture: the conical shape of the Itena posts and its modulus of elasticity, similar to that of dentine, distribute the stresses of mastication evenly along the root, decreasing the occurrence of root fractures. Preserves tooth structure. High mechanical and chemical retentivity. Versatility: Itena posts can be used in narrow, medium and wide root canals. High radiopacity. No corrosion. Saves time and decreases costs. Easy to remove after cementation: the longitudinal arrangement of the fibers favours post removal.

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Use of the dental dam is obligatory. It is essential to retain a sufficient thickness of wall (at least 1 mm) to avoid any perforation or weakening and leave at least 4 to 5 mm of gutta-percha intact to protect the apical seal. Excessive removal of gutta-percha may lead to endodontic failures. The posts are intended for single use. There's re-use may give rise to infection.

The product is intended for practitioners who have received suitable training. If you have any doubts about using this product, consult your local supplier or the manufacturer. Do not sterilize glass fiber posts. Caution: the box is not sterilisable.

CONTRA-INDICATION AND SIDE EFFECTS :

No use on patients with abnormal functional habits (chronic bruxism, ...).

DIRECTIONS FOR USE:

No use on patients presenting a risk of parafunction (bruxism chronic, ...).

CONTRA-INDICATIONS AND EFFETS INDESIRABLES :

No use on patients presenting a risk of parafunction (bruxism chronic, ...).

INSTRUCTIONS FOR USE:

1) Select the correct post according to the diameter and depth of root canal. The choice of the size is ultimately up to the practitioner.

2) Remove part of the root canal filling with pre-warmed instruments to the predetermined depth. 3-5 mm of gutta-percha filling should be left intact to prevent bacterial contamination.

3) Prepare the root canal with the appropriate matching drill for the selected ITENA post (fig 3).

4) Take a periapical radiograph to check the root canal preparation and filling. Excessive removal of gutta-percha may lead to endodontic failures.

Preparation of the post for cementation:

1) Clean the post for cementation.

2) Adjust the length of the post (Itena posts should only be cut with diamond burs at high speed under a flow of water), (fig. 5).

3) Clean Itena posts with alcohol and air-dry.

4) Apply a layer of the selected adhesive according to its instructions (fig 6).

CLASSIFICATION OF ADHESIVES

CONVENTIONAL		SELF-ETCHING	
Three-step	Two-step	One-step	Two-step

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Scellement : The polymerisation reaction of adhesives must be chemical (self-cured) or dual (self- and light-curing).

The reaction of polymerisation of adhesives must be chemical (autopolymerization) or dual (auto-autopolymerization).

The success of the bonding depends on the correct selection of the adhesive.

1) Isolate the area, wash the root canal with air/water spray and dry.

2) Etch the root canal for 15 seconds with 37% phosphoric acid (if you are using a conventional adhesive).

3) Wash the root canal and dry with air and absorbent paper points.

4) Apply the selected adhesive, according to its instructions, into the root canal and onto the post.

5) Mix chemically self-cured or dual-cured restorative composite resin (Dentocore - ITENA) and fill the root canal.

6) Insert the post at the level of the root canal entrance with a diamond disc at high-speed under a flow of water.

7) Reproduce the part of the canal radicular with a spray/air/eau and dry.

8) Place the post in the canal, remove the excess of cement and wait for polymerization.

9) Insert the core with a restorative composite resin (Dentocore - ITENA).

10) Seal the core with a temporary restoration.

11) Prior to making the temporary restoration, isolate the moncone with vaseline.

12) Cut the post at the level of root canal entrance with a diamond disc at high-speed under a flow of water.

13) Make a notch in the top of the post with a round bur at low spe.

14) Place Largo drill running at a low speed onto the notch and the post will disintegrate and be displaced.

15) Realise an encoche in the top of the tenon with a diamond disc at high-speed under a flow of water.

16) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

17) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

18) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

19) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

20) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

21) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

22) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

23) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

24) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

25) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

26) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

27) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

28) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

29) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

30) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

31) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

32) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

33) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

34) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

35) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

36) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

37) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

38) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

39) Place a notch in the top of the tenon with a diamond disc at high-speed under a flow of water.

40) Place a notch in the top of