Ceram-Etch

Hydrofluoric Acid Etching Gel (9%)

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INDICATIONS

• Used for etching of various types of resin bonded porcelain/ceramic restorations (conventional feldspathic porcelain/glass ceramic, lithium disilicate, zirconia reinforced lithium silicate glass-ceramic, fluorapatite glass ceramic, etc.) prior to the bonding and cementation procedure.

MAIN COMPOSITION

• Water, hydrofluoric acid, food grade thickening and coloring agents.

Accessories:

• Dispensing Tips

CONTRAINDICATION

• Not for etching of tooth surface.

OPERATING PROCEDURE

• 1. Never use excessive force on syringe plunger to dispense etch which may lead to over dispensing of the material. If resistance is felt, replace the tip and/or flush excess orally.

• 2. Apply the Ceram-Etch Hydrofluoric Acid Etching Gel onto the entire internal surface of the porcelain restoration, and wait for 60 seconds. One may use the tip to stir the gel in order to refresh the gel contacting the porcelain surface to have a complete etch effect during the waiting time.

• 3. Wash off the Hydrofluoric Acid Etching Gel with a jet of stream water for 10 seconds and dry with contamination-free air.

• 4. Proceed to the application of Silan-It Ceramic Coupling Agent (see the instruction for use of the material).

NOTE:
• Refer to outer package for expiration date. Do not use after expiration date.

• Notification of any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority in their country.

EXTRACTION

• Contactez Hydrofluoric Acid. Avoid contact with skin, eyes and soft tissue. In case of contact with skin or eyes, flush immediately with water. Get medical attention for eyes. Do not take internally.

DISPOSING TIPS

• Dispensing tips are for single use only, to prevent cross-contaminations between patients.

• If there is spill of the Hydrofluoric Acid Etching Gel, pour sufficient amount of base natural chemicals such as Baking Soda, Calcium Hydroxide, etc. to the spilled material to neutralize the acid immediately. Then clean the site with water or wipe it off with proper materials and discard following local environmental regulations.

• Wear acid resistant gloves and eye protection goggles when handling the Hydrofluoric Acid Etching Gel.

• Keep away where children may reach.

STORAGE

• Store Hydrofluoric Acid Etching Gel at 2-25°C and away from direct light. Use the material at normal room temperature.

• Shelf life: 2 years from date of manufacture

WARRANTY

Rizhao HuGe Biomaterials Company, Ltd. warrants this product is free from defects in material and manufacture. HUGE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the products for user's application. If this product is defective within the warranty period, our exclusive remedy and HUGE's sole obligation shall be repair or replacement of the HUGE product.

CONTRIBUTION

• Not to use the product for the mordançage of the surface of the dent.

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