

# Prevent Seal

Self-etching light-cured pit  
& fissure sealant



Self-etching

Adhesion to enamel

Light-curing

Extra-thin tips

Enhance your  
daily practice



# Let's protect pits and fissures



## Key points

### Self-etching

- > No etching or rinsing is required

Adhesion to enamel (20 - 25 MPa)

Contains fluorides

Light-curing to 4.5 mm

### Application with extra-thin tips

- > For application of sealant without material loss

Penetration into grooves and fissures

Extra-thin needle tips  
for necessary application

## Clinical case

Case study by Dr Shalom Melher



Initial view.



1 Clean with a small brush under a stream of water.



2 Dry with care.



3 Apply Prevent Seal.



4 Allow to act for 15 seconds.



5 Light-curing for 20 seconds.



6 Check the occlusion.



Final view.

## Indications

For the preventive sealing of pits and fissures, utilizing the self-etching technique.

Sealing with Prevent Seal gives protection to the susceptible tooth surface by protecting them with impenetrable resin layer, which prevents the penetration of bacteria.

### TIPS AND TRICKS

- > Use air abrasion then rinse the tooth carefully before applying the sealant. No aluminum particles should remain, as they reduce the adhesion of **Prevent Seal**.
- > The prepared surface must be well dried before application of **Prevent Seal**, which is hydrophobic.
- > The **Prevent Seal** is a composite resin, and is compatible with all composites on the market.
- > Contains the camphorquinone photoinitiator. Use a lightcuring lamp which emits between 400 and 500 nm.



### Prevent Seal

1.4 g syringe + 10x 25G needle tips

PVSEAL



Enhance your daily practice

[itena-clinical.com](http://itena-clinical.com)

Central Parc Bat B, 97 Allée de la Louve, 93420 Villepinte - FRANCE - T. + 33 1 45 91 61 40.  
FP.PREVENTSEAL.EN-B-04/2024

Follow us!

