Provitemp

Short-term temporary cement

www.itena-clinical.com

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EN
Non-eugenol, resin based temporary cement containing Fluoride, Chlorhexidine and Potassium Nitrate.

INDICATIONS:
Temporary cementation of crowns, bridges, inlays, onlays and splints.

PROPERTIES:
Provitemp has excellent flow to permit the restoration to be easily and completely seated. • It is strong enough to withstand the stresses of mastication, yet permits easy removal of the restoration. • Provitemp does not contain essential oils so it will not harm temporary crowns or interfere with permanent cementation.

INSTRUCTIONS FOR USE:
1) Complete tooth preparation and fabrication of the temporary restoration before extruding the mixed cement.
2) Dry and isolate the prepared tooth. Be careful not to desiccate or over-dry the preparation.
3) Remove the cap from double barrel automix syringe and set aside. Extrude a very small portion of the unmixed material to insure unobstructed flow from both the base and catalyst components of the respective barrels.
4) Place the mixing tip onto syringe and lock in place with a quarter turn.
5) Extrude a thin layer of the automixed cement onto the cleaned inside surface of the temporary restoration. You have approximately 45 seconds of working time after extruding the mixed cement.
6) Firmly seat the temporary resto- ration in place. Initial gel-set occurs within 1.5 to 2 minutes after seating. Cleaning of excess cement should be initiated immediately after this gel-set phase.
7) Final cement set occurs between 3 and 4 minutes after seating of the temporary restoration.

TECHNICAL DATA:
Working time: 45 s*
Setting time: 3-4 min*
*These indications are obtained at the ideal mixing conditions of 24 °C (75 °F) and 40 % – 60 % relative humidity.

RECOMMENDATIONS/ WARNINGS:
Some individuals may be sensitive to this product. If an allergic reaction occurs, discontinue use and remove from patient. The patient should consult a doctor.

ESSENTIAL INGREDIENTS:
Fluoride • Chlorhexidine • Potassium Nitrate • Methacrylates • Urethan dymethacrylate • Polymerisation activator.

STORAGE:
Keep refrigerated. For dental use only.

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