# Numerys GF

## Glass fiber blocks and disks

### FR

#### DEFINITION

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### EN

#### DEFINITION

NUMERYS GF is a glass fiber and epoxy resin composite for CAD CAM systems. It is available in laboratory disks and unitary blocks.

#### INDICAZIONI

NUMERYS GF is indicated for the preparation of anatomical and aesthetic glass fiber post-and-core prosthesis elements using CAD CAM technologies.

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NUMERYS GF is not indicated for the preparation of other dental restorations (crowns/inlays/onlays/veneers/bridges)

#### ADVANTAGES

- Metal free post and core restorations.
- No core opacification required for higher aesthetics.
- Better mechanical resistance than traditional post and core systems (post and core build-up material)
- Elasticity modulus comparable to dentin's giving the material a flexibility against constraints and minimizing the risk of root fracture.
- Better anatomical fit in root canal promoting a better adhesion between the post-and-core system and the root canal.
- Easy milling of the material using a diamond bur.
- Same root canal and post and core preparation as dental glass fiber posts.
- Same cementation procedure as traditional post and core systems.

#### COMPOSITION

Glass fiber: 75-80%. Epoxy resin: 20-25%

#### APPLICATION

##### MACHINING

1. Scan the preparation of the root directly inside the mouth or the dental wax/ resin model.
2. Select and adjust the material following your usual CAD software system instructions.
3. Position NUMERYS GF in the milling equipment in accordance with the parameters defined by the manufacturer.
4. Proceed to milling of the element using a diamond bur.
5. Remove the post-and-core from the equipment and make the final adjustments using a diamond bur.
6. Dry the post-and-core with air blasting in order to remove all humidity.

#### CLINICAL

Root canal preparation

1. Etch the root canal with phosphoric acid.
  2. Rinse thoroughly and dry with paper points.
  3. Apply your usual adhesive system inside the root canal following the manufacturer's instructions.
  4. Let dry.
- Post-and-core preparation
1. Place the post and core inside the root to check for anatomical fit. If adjustments are required, use a diamond bur.
  2. Clean the post-and-core with alcohol and dry it with air.
  3. Apply a layer of silane on the totality of the post-and-core's surface and let it to dry for 1 min.
  4. Apply your usual adhesive system on the totality of the post-and-core's surface.

NOTE: it is highly recommended to apply silane and adhesive bonding on the post-and-core element before cementation for a perfect adhesion.

5. Cement the post-and-core using a self-adhesive universal resin cement.

NOTE: Do not let the core uncovered after cementation. The impermeabilization of the core part using a resin cement is highly recommended.

Crown preparation

1. Place the crown to check for anatomical fit.
2. Prepare the prosthetic crown following the manufacturer's instructions.
3. Cement the crown using a self-adhesive universal resin cement.

#### PRECAUTION

- Wear a mask and protective glasses during milling and final adjustments of the element.
  - Use a diamond rotary bur for milling and adjustments. Do not use scissors or pliers as it can damage the glass fibers inside the element.
  - Do not use in patients with reported sensitivity to the components inside the formula.
  - This product is exclusively for dental technicians and dentists to use.
  - En ambos casos, la fresa tiene que ser de diamante. No utilice tijeras o alicates, ya que pueden dañar la fibra de vidrio que hay en el interior del elemento.
  - No utilizar en pacientes sensibles a los componentes de la fórmula.
  - NUMERYS GF does not need to be polymerized after processing.
- IMPORTANT: These instructions are based upon clinical and laboratory studies. The success of the procedures depends on the correct diagnosis and careful technique.

#### STORAGE

Maintain at room temperature and protected from light.

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