# PREVENT SEAL

FR

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EN

Self-etching light cured, pit & fissure sealant with fluoride release.

#### INDICATIONS:

For the preventive sealing of pits and fissures, utilizing the self-etching technique. The pits and fissures of teeth are particularly at risk of dental tooth decay since they favor the colonization of microorganisms. Sealing with Prevent Seal gives protection to the susceptible tooth surface by protecting them with impenetrable resin layer, which prevents the penetration of bacteria.

#### PROPERTIES:

Single component material. Self-Etching: NO ETCH / NO RINSE / NO DRY. Single Step. Easy & precise application by using direct syringe technique. Light-cured. Low viscosity, excellent flow properties, optimal wetting of fissures. Fluoride release.

#### INSTRUCTIONS FOR USE:

- 1) Clean the enamel surface thoroughly with either cleansing paste containing neither oil nor fluoride, or with pumice and brush.
- 2) Rinse thoroughly and air dry. Avoid contamination of dried surface with saliva.
- 3) Mount and secure lock a new dispensing tip on the Prevent Seal syringe and carefully syringe an adequate amount of Prevent Seal into the fissures. Let the sealant flow until the entire fissure is covered. You can use a disposable brush, micro-brush or the dispensing tip itself to spread the sealant and eliminate air bubbles.
- 4) Allow Prevent Seal to sit on the enamel surface for 15 sec.
- 5) Light cure for 20 sec.
- 6) Check occlusion.

#### CONTRAINDICATIONS:

The sealing of wells and grooves with Prevent Seal are not effective if the following occur: when wells and grooves have already been attacked by caries. If the affected working area is not kept dry during the application and photopolymerisation of Prevent Seal. If the patient has an allergy to any or part of the ingredients of Prevent Seal.

#### CAUTION!

Prevent Seal contains polymerizable monomers, which may cause allergic contact dermatitis in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if known allergy to methacrylate resin exists, discontinue use.

Recommended Wavelength: 465 nm. Intensity of the light source recommended: 600mW/cm<sup>2</sup>.

#### ESSENTIAL INGREDIENTS:

Urethandimethacrylate oligomer - Bis-GMA, Triethyleneglycoldimethacrylate - 2-Hydroxyethylmethacrylate - Self-etching agent - White pigment - Glass fillers - Functional monomers - Urethan Dimethacrylate - Triethyleneglycoldimethacrylate - 4-Methacryloxyethyltrimellitic anhydride - Photoinitiators - Fumed Silica.

#### COMPONENTES FUNDAMENTALES:

Oligómero de dimetacrilato de uretano - Bis-GMA - Trietilenoglicoldimetacrilato - 2-hidroxietilmetacrilato - Agente de auto-aplicação - Pigmento branco - Rellenadores de fibra de vidro - Monómeros funcionais - Dimetacrilato de uretano - Trietilenglicol dimetacrilato - 4-Metacriloxietilmelítico anhídrido - Fotoiniciadores - Silica pirogénica.

#### STORAGE:

Do not store above 24°C (75°F). Keep away from direct sunlight.

For optimum freshness, keep refrigerated. Do not freeze. If refrigerated, allow the syringe to reach room temperature. Do not use after expiration date. Shelf life is 2 years after manufacturing date.

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