# Reflectys Flow

## Flowable universal composite

### FR

#### Composite fluide micro-hybride nano-renforcé photopolymérisable fluoré

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### EN

#### Flowable and fluorized nano-reinforced micro-hybrid light-cured composite

**DESCRIPTION :**  
REFLECTYS Flow is a light-cured, radiopaque flowable restorative composite resin.

**COMPOSITION :**  
Bis-GMA - TEGDMA - UDMA - Barium glass - Silica - The content of inorganic filler is ca.60% by weight.

**INDICATIONS:**

- Class V restorations (cervical caries, root erosion, wedge shaped defects)
- Anterior restorations (Class III, IV)
- Small posterior restorations
- Extended fissures sealing in molars and premolars
- Repair of composite/ceramic veneers
- Blocking out of undercuts

**CONTRA-INDICATIONS:**

REFLECTYS Flow should not be used with patients who develop hypersensitivity to it or have a history of hypersensitivity to methacrylate monomer.

**RECOMMENDED PROTOCOL:**

1. Shade selection: Clean tooth with pumice powder/paste and water to remove surface stains or extraneous oxide. Prior to isolation of tooth, select the appropriate shade.
2. Cavity preparation: Remove all amalgam or other base materials that interfere with operations.
3. Pulp protection: In deep cavities cover the dentin close to the pulp with a minimum amount of calcium hydroxide liner leaving the rest of cavity surface free bonding. Glass ionomer or other eugenol-free materials may be used, if wished.
4. Enamel and dentin treatment: Follow the manufacturer's instructions regarding etching, priming, adhesive application and curing.
5. Dispensing: After removing cap, insert the disposable tip onto the syringe. Before applying resin make sure that needle of tip is not plugged. If plugged, remove the tip and express a small amount of resin directly from the luer of syringe. Remove any visible plug, if present from the syringe. Replace new disposable tip and again express resin. REFLECTYS Flow may be extruded onto a dispensing container and applied with a brush or other instruments.
6. Placement: Place resin in increments in 2.0mm levels or less. To prevent oozing, pull back on syringe plunger after dispensing procedure.
7. Curing: Expose each area of restoration surface to a high intensity visible light source. Hold the light guide tip as close to the restorative as possible during exposure. Light cure for 20 seconds with standard light curing unit. (If light curing unit output is below 400mW/cm², as measured by a curing radiometer, more time may be needed.)
8. Finishing: finish and polish using conventional techniques.

**STORAGE:**

1. Do not store at elevated temperature or intense light.
2. Keep in a cool place.
3. Do not store material in proximity to eugenol containing product.
4. After using REFLECTYS Flow, it should be capped immediately.

**PRECAUTIONS:**

1. For dental use only.
2. It should be used by a dentist and dental professional. Observe the procedure and the precautions.
3. Keep it away from a reach of non-dental user including children, senior citizens and others.
4. If skin sensitization or allergy to acrylate occurs, stop using resin.
5. Exposure to uncured resin should be avoided. If contacts skin, wash immediately with soap and water.
6. If accidental contact with eyes or prolonged contact with oral tissues occurs, flush immediately with a large amount of water.
7. If refrigerated, the material must be allowed to reach room temperature before use.
8. Use of protective gloves and a no-touch technique is recommended.
9. Protective device for the respiration tract is highly recommended before using.
10. In case of non-polymerization, check the light intensity of the light curing unit and then contact the nearest dealer(s) or manufacturer directly regarding the same quality matters.
11. Do not use after expiry date.

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