FR

INDICATION

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
EN

INDICATION

Posts and keys are designed to support a temporary or permanent reconstitution using a direct or indirect method, following root treatment. Drill bits are designed to prepare the posts or key holes.

INSTRUCTIONS FOR USE:

- Choose the suitable Denticloic pin by superimposing the calibration film on the periapical X-ray.
- Prepare the root canal by using:
 1. The Gates or Largo drill bits: remove the gutta-percha (A)
 2. The Denticloic «Pilot» drill bits: help widen the canal.
 3. The Denticloic cylindrical-conical or conical drill bits: prepare the holding place for the root post according to the practitioner's choice. (B)
- Check that the root is whole before inserting the post.
- To make things safer, use drill bits of increasing diameter in succession to widen the canal progressively to the desired diameter.
- Ensure that sufficient thickness of dentine is retained. (C)
- Insert the metallic pin for taking the impression. Send to the laboratory with the pin in the impression, and attach the corresponding burnout post.
- In cases where the roots diverge, attach the corresponding Denticloic posts).
- Insert a metallic pin of the same diameter into the canal and produce a temporary crown. Place the temporary cement on the edges of the crown only, not on the post.
- Cast the master model in the laboratory.
- Insert the burnout post into the master model and produce the model. This model is cast using the conventional method.
- Then test the prosthesis and apply sealing for definitive restoration. (D)

PRECAUTIONS FOR USE:

- Do not use the burnout posts for taking the impression. There is a 2/100 difference between metallic pins and burnout posts; burnout posts are thinner in order to provide for the thickness of the sealing cement and to prevent excessive friction between the metal and the tooth.
- Pay close attention to the lengths and diameters of the pins.
- Denticloic stainless steel and titanium pins must be sterilised before use.
- Denticloic drill bits must be sterilised before the first use and after each use.
- Remove the colour code before inserting the post.
- The Denticloic burnout posts and the principal box cannot be sterilised. Only the Denticloic module may be sterilised.
- The sterilisable Denticloic module allows drill bits and posts to be sterilised before each use. It is advisable to place two drill bits and two posts only in the sterilisable module.
- Use of the dam is essential.
- Denticloic drill bits must be replaced after every 15 uses.
- Cleaning and sterilisation: All metallic parts may be cleaned using ultrasound and sterilised in an autoclave (in our sterilisable kit). The posts and instruments must be cleaned prior to use, by immersion in a disinfectant and decontaminant solution, and then sterilised (134°C for 18 minutes, continuous temperature), using a validated process available on the Itena Clinical website: <http://www.itena-clinical.com/en/>
- Denticloic drill bits are designed to be used with standard contra-angles at slow speed, without exceeding 10,000 rpm.

HYGIENE: CLEANING/STERILISATION PROTOCOL

Recommendations

- Sterilisation of the posts before first use or before each subsequent use is the responsibility of the user. Itena Clinical cannot be held liable for use of dirty and/or damaged instruments by the user.
- Do not reuse single-use devices (Itena posts and screw-posts).
- Drill bits are designed for a maximum use of 15 cycles.
- Clean the instruments as soon as possible after use.
- Do not use detergent solutions with an acid, abrasive or chlorine base.

Manual cleaning

- Clean the instruments with a suitable nylon brush in order to remove dentine residue. No traces of dirt may remain on the posts.
- Immerse the posts completely for 15 minutes in a disinfectant solution suitable for medical devices*.
- Ensure that no dirt is left on the posts. If necessary, repeat the cleaning and disinfection cycle.
- Rinse the posts thoroughly with demineralised or distilled water to prevent corrosion and alteration of shape during the sterilisation process. Dry with compressed air for 2 minutes. Check visually for cleanliness and discard any instruments that may be distorted.

Automated cleaning

- Use a cleaning and disinfection device, sorting the instruments according to type and category of material in order to avoid electrolysis and corrosion.
- Position the posts so that they are completely immersed, and launch the cleaning cycle according to the following indications:

Stages	Type of product and concentration	Time / Temperature
Pre-wash	/	/
Washing	Disinfectant solution suitable for medical devices*	15 minutes at 55°C
Rinsing	Demineralised water	2 minutes at ambient T° (15-25°C)
Drying	/	According to autoclave parameters

- At the end of the cycle, ensure there is no dirt left on the posts. Repeat the cycle if necessary.

STERILISATION

- Place the posts in suitable sterilisation sachets (compliant with EN ISO 11607).
- Place the sachet in the autoclave and then perform a sterilisation cycle: 134°C at 2.2 bars for 18 minutes.
- Keep the posts away from any source of contamination.

Warning!

GENERAL RISKS:

This product is designed for suitably trained professional persons. If you have any doubts regarding use of the product, contact us. Because of the risk of the post or an accessory falling into the mouth, use a dam. The choice of size, shape and material of the post is the responsibility of the practitioner. Risk of breakage of burnout posts: never leave these inside the mouth. The compatibility and safety of posts have not been assessed in a magnetic resonance environment. Ask the patient to make systematic provision when magnetic resonance examinations are performed.

RISKS OF INFECTION:

To avoid the risk of infection, follow the advice relating to sterilisation. Reuse of posts can carry a risk of infection. Used products must be disposed of as potentially infectious medical waste.

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