99 mm

99 mm

99 mm

ITENA®

DENTO CLIC

GLASS FIBER

Dentoclic Glass Fibre Post

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
GLASS FIBER

Dentoclic Glass Fibre Post

Glass fibre posts for the reconstruction of endodontically treated teeth.

COMPOSITION (% weight):

• Posts Glass fibre 80% • Epoxy resin 20%.

INDICATIONS:

Prosthetic reconstruction of a tooth following root treatment.

PROPERTIES:

• Aesthetics. • Reduces the risk of root fracture: the conical shape of Itena posts and its modulus of elasticity, similar to that of dentine, distributes the stresses of mastication evenly along the root, decreasing the occurrence of root fractures. • Preserves tooth structure. • High mechanical and chemical retentivity. • Versatility: Itena posts can be used in narrow, medium and wide root canals. • High radiopacity. • No corrosion. • Saves time and decreases costs. • Easy to remove after cementation: the longitudinal arrangement of the fibres favours post removal. • The flexural strength of Itena posts is 857 MPa, according to ISO Standard 10477.

WARNING:

• Use of the dental dam is obligatory. • It is essential to retain a sufficient thickness of wall (at least 1 mm) to avoid any perforation or weakening and leave at least 4 to 5 mm of gutta-percha intact to protect the apical seal. • Excessive removal of gutta-percha may lead to endodontic failures. • The posts are intended for single use. There re-use may give rise to infection and/or cross-contamination. • This product is intended for practitioners who have received suitable training. • If you have any doubts about using this product, consult your local supplier or the manufacturer.

Caution: the box is not sterilisable.

CONTRA-INDICATION AND SIDE EFFECTS:

• Do not use on patients with abnormal functional habits (chronic bruxomania, ...).

DIRECTIONS FOR USE:

Root canal preparation: 1) Select the correct post according to the diameter and depth of root canal. Use the ruler and the periapical radiograph of the tooth (fig 1) - The choice of the size is ultimately up to the practitioner. 2) Remove part of the root canal filling with pre-warmed instruments to the predetermined depth. 3- 5 mm of gutta-percha filling should be left intact to prevent bacterial contamination (fig 2). 3) Prepare the root canal with the appropriate matching drill for the selected Itena post (fig 3). Drills must be replaced every 40 to 60 uses. 4) Take a periapical radiograph to check the root canal preparation and filling. Excessive removal of gutta-percha may lead to endodontic failures.

Preparation of Itena post for cementation: 1) Check insertion and position of the post in root canal (fig 4). 2) Adjust the length of the post (Itena posts should only be cut with diamond burs at high speed under a flow of water (fig 5)). 3) Clean Itena posts with alcohol and air-dry. 4) Apply a silane agent (e.g. Silanea - Itena) and let it dry. 5) Apply a layer of the selected adhesive according to its instructions (fig 6).

Cementation: The polymerisation reaction of adhesives must be chemical (self-cured) or dual (self- and light-cured). Successful post cementation is highly dependent on adhesive selection.

CLASSIFICATION OF ADHESIVES			
CONVENTIONAL		SELF-ETCHING	
Three-step	Two-step	One-step	Two-step
Phosphoric acid	Phosphoric acid	Acid	Acid
Primer	Primer	Primer	Primer
Adhésif	Adhesive	Adhesive	Adhesive

1) Isolate the area, wash the root canal with air/water spray and dry. 2) Etch the root canal for 15 seconds with 37% phosphoric acid (if you are using a conventional adhesive). 3) Wash the root canal and dry with air and absorbent paper points. 4) Apply the selected adhesive, according to its instructions, into the root canal and onto the post. 5) Mix chemically self-cured or dual-cured restorative composite resin (e.g. Dentocore - Itena) and fill the root canal. 6) Insert the post in the canal, remove the excess cement and pause to allow polymerisation. If the restorative composite resin is dual-cured, light-cure it at this point (fig 7).

Building the core (coronal part): 1) Build the core with a restorative composite resin (e.g. Dentocore - Itena). 2) Reproduce the coronal portion (fig 8). 3) Prior to making the temporary restoration, isolate the core with petroleum jelly.

Technique for the eventual removal of post (in case of necessity): 1) Cut the post at the level of root canal entrance with a diamond disc at high-speed under a flow of water. 2) Make a notch in the top of the post with a round bur at low speed. 3) Place a Largo drill running at a low speed onto the notch and the post will disintegrate and be displaced.

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