DENTOCORE *body*

CORE BUILD-UP

DentoCore Body Core Build-up

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CORE BUILD-UP

DentoCore Body Core Build-up

Nanofill composite core build-up material. Dual-cure.

INDICATIONS

DentoCore Body is a dual-cure, fluoride releasing, resin composite material designed for the fabrication of core build-ups. DentoCore Body may be used with any Bis-GMA compatible bonding agent. DENTOCORE and DENTOCORE Body having the same properties and indications, the choice of the material depends on dentist habits.

PROPERTIES

• Easy to use. • Provides a consistent, homogeneous mix. • Dual curing: material can be light cured, additional chemical cure assures that the material can be quickly, bulk filled. • Dense void-free core: cuts like dentin. • Fluoride release: strengthens natural tooth structure and helps prevent secondary caries. • Radiopaque. • Superior compressive strength for long-lasting restorations. • Enhanced mechanical properties of conventional Bis-GMA composite by utilization of Intercalated and Exfoliated nanocomposite technology.

INSTRUCTIONS FOR USE

1) Prepare the tooth according to clinical requirements. 2) Isolate the tooth from oral fluids. A rubber dam is recommended. 3) Acid etch the entire preparation with a 32%-40% phosphoric acid for 15 seconds or apply a self etching adhesive and follow instructions for use of the manufacturer. 4) Rinse off etchant with a water spray for at least 15 seconds and apply bonding agent of choice. 5) Place into base of prepared tooth an adequate amount of DentoCore Body to complete build-up. 6) Do not disturb DentoCore Body during the gelation stage. Light-cure for at least 40 seconds. Check for adequate hardness. Chemical-cure intraoral setting time is approximately 4.5 minutes from start of mix/application. Set time will be longer outside of mouth.

TECHNICAL DATA

Working Time at 23°C	1 : 30 minute
Setting Time at 37°C	4 : 30 minutes
Compressive strength	250 MPa
Diametral tensile strength	40 MPa
Flexural strength	200 MPa
Water sorption	14 µg/mm ²
Solubility	2 µg/mm ²
Linear shrinkage	1,2%
Radiopacity	400 %Al
Compatible with	halogen light, plasma arc light and LED

CAUTION!

DentoCore Body contains polymerisable monomers which may cause allergic contact dermatitis in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if known allergy to methacrylate resin exists, discontinue use.

PRECAUTIONS IN USE

• Do not use in presence of eugenol based materials; complete setting of DentoCore Body may not take place. • A calcium hydroxide or glass ionomer liner may be used to isolate DentoCore Body from eugenol sources. • If refrigerated, allow to return to room temperature before use. Refrigeration will increase viscosity, making application more difficult, and will retard setting time. • Single use tips.

ESSENTIAL INGREDIENTS

Bis-GMA, TEGDMA, Silica, Barium aluminoborosilicate, photo-initiator, cross-linker

STORAGE:

• Do not store above 24°C (75°F). • For optimum freshness, keep refrigerated. • Do not freeze. • Do not use after expiration date. • See expiration date on the label.

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