**FR**

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**EN**

**DESCRIPTION:**  
Dentocrown HD cartridge into a dispenser with a suitable piston. Due to technical reasons, slight differences in the filling level between the cartridge chambers can occur.

**COMPOSITION:**  
Bis-GMA - Bis-EMA - TEGDMA - UDMA - micro/nano silica - others

**INDICATIONS:**  
• Fabrication of temporary crown, bridges, inlays, onlays, partial crown and veneers.

**APPLICATION:**  
Insert Dentocrown HD cartridge into a dispenser with a suitable piston. Due to technical reasons, slight differences in the filling level between the cartridge chambers can occur.

**INDICATIONS FOR THE APPLICATION:**  
First, fill the retentive parts in the preparation, do not use cements containing eugenol.

1. Prepare an impression prior to extraction of a tooth in cases where no prefabricated crowns or similar restoratives are used.

2. Dry the teeth and lightly lubricate the area by applying a lubricant such as Vaseline or similar materials.

3. Dispense Dentocrown HD directly into the impression or matrix. To prevent bubbles, it is important to keep the tip of the mixing cannula immersed in the material. Do not overfill.

4. Gently insert into the mouth and firmly hold the prepared tooth (within its impression matrix) for up to 45 seconds.

5. Approximately 1:30 to 1:45 minutes from start of mix, the provisional crown or bridge can easily be removed together with the impression from the prepared teeth. The setting reaction has to be monitored intraorally, as the temperature in the mouth has a significant influence on the setting time and removal of the provisional is only possible during gel state.

6. The crown and bridge may be contoured and polished after 1 additional minute. The removal of the soft, sticky inhibition layer on the surface of Dentocrown HD (e.g. by wiping with ethyl alcohol) facilitates the shaping. Completed crown can be finished using toothed grinders, etc.

7. Carefully trim any excess material and roughen area around gingiva. Cement the temporary restoration using an eugenol-free temporary cement.

**CAUTION:**  
• This material is designed for sale and use by dental care professionals only.

• This material contains methacrylic monomer. Some patients can have an allergic reaction to the monomer. Do not use this material for patients allergic to or hypersensitive to methacrylic monomers.

• If this material causes an allergic reaction or oversensitivity, discontinue its use immediately.

• Avoid contact with eyes, mucosal membrane, skin or clothing.

• If this material comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.

• Use a rubber dam to protect the throat.

• Do not use this material after the indicated expiration date on the cartridge.

• If refrigerated storage is used, allow product to reach room temperature prior to use.

• Do not use when pulp exposure has occurred and the situation warrants a direct or indirect pulp capping procedure.

• Single use tips.

**STORAGE:**  
1. Do not expose restorative materials to elevated temperature or intense light.

2. Store this material at a temperature between 2-27°C

3. Ambient temperature routinely higher than 27°C/80°F may reduce shelf life.

4. Shelf life at room temperature is 24 months.

5. Use over cartridge for expiration date.

6. Do not store materials in proximity to products containing Eugenol.

**WARRANTY:**  
Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the suitability of the materials delivered by us for the intended purposes of application. As the application of our preparation is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

**CUSTOMER INFORMATION:**  
No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

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