99 mm

ITENA

# DENTOTEMP

LONG-TERM TEMPORARY CEMENT

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# DENTOTEMP

LONG-TERM TEMPORARY CEMENT

## DentoTemp Long-Term Temporary Cement

Temporary sealing cement with urethane-acrylic polymer base, free of eugenol. Specially formulated for sealing of prosthetic elements on implants.

### INDICATIONS:

- Temporary sealing on natural, non-reconstituted teeth (composite, core).

FOR RECONSTITUTED TEETH, USE Provitemp.

- Temporary sealing for cosmetic elements such as veneers or onlays.
- Permanent sealing for crowns and bridges on implants with possibility of removing restorations.

### PROPERTIES:

- DentoTemp is a temporary sealing cement free of eugenol, zinc oxide and essential oils, which allows ease of cleaning and prevents interference during the final and permanent adhesion of the prosthetic part. Polymerisation is carried out in two stages to allow faster, simpler cleaning, without interfering with work time or compromising performance.

- DentoTemp is ideal for permanent sealing of crowns on implants. It provides an easy seal and a high level of retention, while keeping open the possibility of removing restorations, if necessary.

- Low levels of solubility in mouth fluids, and rapid setting, mean a high rate of retention.
- Excellent stability is provided in the mouth, with a very good mechanical resistance.

- No taste and no smell mean more comfort for the patient.

### INSTRUCTIONS:

Mixing of DentoTemp manually:

- 1) Apply equivalent measures of base and catalyst to the mixing block.
- 2) Mix them thoroughly with a spatula for 20-30 seconds.

Automix :

Purge the syringe before installing the mixing tip on the syringe and purge the syringe again.

### Sealing of the temporary fixed prostheses:

Before sealing, check the adjustment, occlusion and proximal contact points, and carry out all the finishing procedures. 1) Dry and isolate the prepared surfaces. Absolute drying is not essential, but contamination with saliva must be avoided. 2) Ensure that the intrados of the temporary crown is clean and dry. Apply the cement mixture directly into the temporary crown. Working time: 45-60 seconds. 3) Put the temporary crown into position immediately, hold it in place or have the patient bite gently on a cotton roll until initial setting occurs (2-2½ minutes). Initial setting is characterised by a distinct rubbery phase during which the excess cement can be easily removed using a probe. 4) Take care to pass the probe through the inter-proximal embrasure to remove excess cement from these areas, before final setting time. 5) Setting time within the mouth is about 4-5 minutes, depending on conditions.

### Permanent sealing of permanent crowns on implants:

The sealing procedure for a permanent crown on a metallic or metal-ceramic implant on an implant abutment is identical to the sealing of a temporary crown. It is not advisable to use a titanium etching agent on the pillars before sealing occurs. If the crown has to be removed, a gentle tap with a hard ring should be sufficient.

### PRECAUTIONS:

Although this is rare, some patients may be sensitive to this product. If an allergic reaction occurs, stop using the product, remove it from the patient's mouth and refer the patient to a doctor. Single use tips.

### COMPOSITION:

- Bis-GMA
- HEMA
- TEGDMA
- Silicon - Other Substances

### STORAGE:

- For optimum freshness, keep refrigerated.\* (temperature below 24°C)

\* If the syringe is refrigerated, it must reach room temperature before being used.

- For dental use only.

99 mm

ITENA

# DENTOTEMP

LONG-TERM TEMPORARY CEMENT

## DentoTemp Long-Term Temporary Cement

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