MTA Bioseal

Bioceramic root canal sealer

FR

DESCRIPTION DU DISPOSITIF

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EN

DEVICE DESCRIPTION

MTA BIOSEAL is a bioceramic root canal sealer based on mineral trioxide aggregate (MTA). It comes in a package containing a 4 g syringe of MTA sealer, 10 mixing tips, and 1 mixing pad.

COMPOSITION

The product consists of two pastes:

- Base Paste: Calcium Tungstate, Salicylate Resin, Natural Resin, Silicon Dioxide and Titanium Dioxide.
- Catalyst Paste: Sulfonamide Resin, Mineral Trioxide Aggregate, Silicon Oxide, Natural Resin and Titanium Dioxide.

INTENDED PURPOSE

Intended use

MTA BIOSEAL is a bioceramic root canal sealer based on MTA (Mineral Trioxide Aggregate) intended for the permanent sealing of root canals and may be used in combination with root canal obturation materials.

Target population

Patients who present permanent teeth with closed apex, diagnosed with pulp and periapical disease, and requiring endodontic procedures like permanent sealing of root canals.

Intended users

MTA BIOSEAL is intended to be used by dental healthcare professionals trained to perform endodontic procedures.

CONTRAINDICATIONS

MTA BIOSEAL contains components that may cause reactions in sensitive individuals. Do not use on patients allergic to resins or other components of the product.

WARNINGS

- Avoid contact with eyes or skin. In case of contact, rinse immediately with water.
- Avoid contact with oral mucosa. In case of contact, rinse with water and prevent swallowing of product. In case any sensitivity persists, seek medical attention promptly
- Any use other than stated indications is under the responsibility of the dentist.

PRECAUTIONS

- When using this product, it is recommended that patients and professionals wear suitable protective clothing, eye protection and gloves.
- The mixing tips and mixing pad are intended for single use only. Disinfection and/or sterilization may affect the characteristics of the product, resulting in additional risk of physical injury or infection to the patient.
- If the syringe becomes contaminated with saliva or blood during application, discard the syringe and do not use it on other patients.
- Please check the integrity of the device before using this device (packaging and syringe).
- Do not use it if the product is damaged.
- Ensure that the caps of the base paste and catalyst paste are not switched, as this may cause the product to harden in the tube.
- The safety and efficacy of this product have not been established in pregnant, lactating women or for children and adolescents.
- The device must only be used by trained and experienced dentists.
- The working length in the root canal should not exceed 2 mm from the radiological apex to avoid apical extrusion.
- The evaluation of the allergic background of a patient is the healthcare professional's responsibility.
- If the MTA BIOSEAL notice is downloaded and saved locally (e.g. on a phone or tablet), it is the dentist's responsibility to ensure that they are the latest version. The most recent version is available on the ITENA Clinical website

OPERATING PROCEDURE

1. Root canal preparation: It is essential to perform endodontic treatment by adhering to a precisely determined working length and in accordance with established protocols. Prior to insertion of MTA-BIOSEAL, the root canal should be prepared, cleaned and dried using paper points, according to the selected endodontic technique. The device requires an isolated dry field; therefore, moisture found in the canal or in the ambient humidity is sufficient to initiate MTA-BIOSEAL hardening;
2. Mixing: The dual syringe ensures equal mixing in a 1:1 ratio. Press plunger to extrude material directly onto a glass slab, mixing pad or into the tooth canal. The mixture should be used immediately after mixing;
3. Insertion: Use MTA-BIOSEAL with gutta-percha coating them with a thin layer of the sealer and filling the root canal according to the selected technique.

WARNING: The self-mixing tip must be discarded after use in accordance with local regulations.

Working time: 23 minutes; Setting time: until 150 minutes

4. Take a radiograph to verify the correct filling of the canal.

REMOVAL

Removal of the root canal filling: Use the conventional techniques for the removal of gutta-percha fillings.

DISPOSAL

The product must be disposed of in accordance with local regulations.

UNDESIRABLE SIDE EFFECTS AND RESIDUAL RISKS

Allergic reaction and temporary irritation of surrounding tissues.

FOLLOW-UP

A post-operative radiograph must be taken to verify that the endodontic obturation has been performed perfectly.

EXPECTED LIFETIME

At least 5 years.

SPECIAL STORAGE CONDITIONS

- Keep in dry and cool place, away from humidity and sources of heat. Do not store the product above 30°C.
- Do not store in the refrigerator
- Keep the product tightly capped.
- Do not use after expiry date.

CLINICAL BENEFITS AND PERFORMANCE

- Adequate radiopacity.
- Suitable flow, allowing fill the accessory and lateral canals.
- Highly alkaline, difficulty of bacterial growth.
- Radiopacity adequate, allowing radiographic visualization.
- Ability to aid bone recovery and cementum formation.

REPORTING OF INCIDENTS

If an accident occurred, report it to the manufacturer and the relevant authority of the concerned country.

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