FR

INTRODUCTION :

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EN

INTRODUCTION:

MTA – Mineral Trioxide Aggregate – has established its use in Endodontics due to the clinical results proven through numerous scientific studies. However, due to the physical properties inherent to the product powder, the “sandy” consistency hinders manipulation, reproducibility after mixing and the transport of material to the site of repair. To improve these characteristics, ITENA has developed a new formulation called MTA BIOREP. This new formula preserves all the chemical and biological properties of the original MTA, which guarantees success of the treatment, but changes its physical handling properties. This results in a product with a greater plasticity, simplifying manipulation and the insertion into the dental cavity.

DEFINITION:

Bioceramic Reparative cement: MTA BIOREP is an endodontic bioceramic reparative cement with high plasticity. It is composed of mineral oxides in the form of thin hydrophilic particles.

INDICATIONS:

- Treatment of root canal and furcation perforations, caused by carious or iatrogenic lesions
- Treatment of root perforations, caused by internal resorption
- Pulp capping
- Pulpotomy (removal of the affected crown portion of the pulp, preserving the vitality and function of the remaining radicular portion)
- Apexogenesis (induction of root development in vital teeth with an inflamed coronal pulp)
- Apexification (induction of a mineralized barrier formation at the root tip of young permanent teeth with incomplete root development and a necrotic pulp)
- Periapical surgery with reverse filling

ADVANTAGES:

- The particles’ size after hydration in this new formula allows an easy manipulation and insertion into the dental cavity.

- The addition of a new radiopacifier, Calcium Tungstate (CaOW4), does not cause root or crown staining and discoloration.

- Low solubility

- Excellent marginal sealing ability through a controlled setting expansion, which hinders the migration of microorganisms and fluids into the root canal.

- Excellent biological sealing of root perforation sites due to peri-radicular cement formation.

- The induction of a dentine shield formation when used on pulp exposures in cases of pulp capping.

- Finally, MTA BIOREP can be used in a wet environment without an alteration of its properties

COMPOSITION :

Powder :
Tricalcium silicate 3CaO.SiO2
Dicalcium silicate 2CaO.SiO2
Tricalcium aluminate 3CaO.Al2O3
Calcium oxide CaO
Calcium tungstate CaWO4

Liquid:

- Water
- Plasticiser

PROPERTIES:

Setting time: MTA BIOREP solidifies when kept in a wet environment after spatulation and compaction. The initial setting time is approximately 15 minutes

Radiopacity: Similar to that of gutta-percha and more radiopaque than dentine or bone. Absence of dental discoloration is due to the new radiopacifier used (CaWO4)

APPLICATION:

Warning: When opening the capsule, please ensure that it stays in a vertical position to avoid any material loss. (See illustrated protocol)

Automatic mix:

- Sterilize all material used for MTA BIOREP placement, insertion and compaction.
- Delicately open the single-dose capsule like illustrated in the protocol.
- Pour 4 drops of the liquid mix directly inside the capsule.
- Close the capsule and place it inside a standard mixing device (around 4200 rpm) for 20-30 seconds. (The final mix should be homogeneous and very workable/malleable)
- Place MTA BIOREP on the desired area with an instrument of choice and condense it with an instrument suited for this purpose.

Manual mix:

- Sterilize all material used for MTA BIOREP placement, insertion and compaction.
- Delicately open the single-dose capsule like illustrated in the protocol.
- Dispense the content of the single-dose capsule on a mixing pad.
- Pour 4 drops of the liquid mix on the mixing pad.
- Spatulate for 1 minute until the powder and liquid are completely homogenized. (The final mix should be homogeneous and very workable/malleable)
- Place MTA BIOREP on the desired area with an instrument of choice and condense it with an instrument suited for this purpose.

PRECAUTIONS:

Use eye protection, a mask and gloves when handling MTA BIOREP. In the case of contact with eye or skin, wash abundantly with water. The product is very sensitive to surrounding humidity, please make sure to open MTA BIOREP capsules only right before use. Do not use MTA BIOREP to fill a canal, because its viscosity and flow are unsuitable for this application. Do not use MTA BIOREP on a tooth surface that is in contact with the gingival sulcus because it will lead to the product’s disintegration. Only use MTA BIOREP after having controlled the acute phase of the endodontic pathology. MTA BIOREP’s solidification is altered by acidic pH of endodontics lesions and surrounding areas. When applying MTA BIOREP, avoid letting it overflow to farthest regions beyond the dental cavity. Even if the excess is reabsorbed, it can slow the healing process.

Warning: The information provided in this manual is based on laboratory and clinical studies. The successful use of MTA BIOREP depends on the correct diagnosis and associated operative technique, the condition of the treated tooth and the general health of the patient. This product should be used according to the instructions described in this manual. The manufacturer is not responsible for the errors or damages caused by an incorrect use of this product or by its use in a non-compliant situation.

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