NUMERYS | HC

 **iTENNA[®]**

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Hybrid ceramic blocks and disks, compatible with all dental CAD/CAM systems.

DESCRIPTION

NUMERYS HC, in the form of hybrid ceramic blocks and disks, is a dental restorative material used to make inlays, onlays, veneers and crowns.

COMPOSITION

NUMERYS HC is a hybrid ceramic material with excellent aesthetic results. It is composed of ca. 80% of mineral filler with particle size ranging between 10nm and 3 μ m.

PROPERTIES

NUMERYS HC has shock absorption properties that ensure better resistance to masticatory forces. The material is also radiopaque so that it can be viewed on x-rays once placed in the mouth.

INDICATIONS


NUMERYS HC is indicated for the indirect fabrication of metal-free prosthodontic components:

Inlays, Onlays, Veneers, Crowns, Implant-supported crowns


CONTRAINDICATIONS

Do not use this material in patients with known methacrylate allergy.


RECOMMENDED PROCEDURE


CROWN

ONLAY

INLAY

VENEERS

1. Prepare the tooth according to the design recommendations detailed below.

Designing restorations fabricated with Numerys HC

- Prostheses milled with Numerys HC follow the standard procedures for CAD/CAM ceramic restorations.
- Round the inner and marginal edges.
- It is advisable to round the shoulder and/or the inner edges of your preparation.
- Each value shown in the diagram shows the recommended minimum thickness for any restoration fabricated with Numerys HC hybrid ceramic.

MILLING

1. Follow the instructions of the milling unit being used.
2. Select the appropriate shade and size of block / disk to match the wanted restoration.
3. Select the appropriate milling program for hybrid ceramic materials.
4. After milling, check the integrity of the piece (discolouration, fissures, cracks).
5. NUMERYS HC does not require a post-milling firing phase.

BONDING

1. Test the milled prosthesis in the mouth so that, if necessary, it can be adjusted with polishers and diamond burs under spray.
2. Before bonding, apply a primer for ceramics (C-Ram Booster) or etch the surface with hydrofluoric acid (Ceram-Etch), then apply a universal silane (Silan-it) on the inner section of the intrados restoration.
3. Then apply a permanent, self-adhesive, dual-cure composite cement (Total C-RAM) and place the definitive restoration on the tooth/teeth.
4. Light-cure for 2 seconds to remove excess cement, then light-cure each surface for 20 seconds.
5. Finish by polishing, if necessary.
6. If necessary, adjust the shade using a light-cured characterisation material

STORAGE

1. Avoid direct exposure to sunlight.
2. Store at room temperature.

PRECAUTIONS

1. NUMERYS HC should be used by a dental professional.
2. Keep out of reach of people not trained to use this product.
3. Use of protective gloves and glasses is recommended.
4. Do not use after the expiry date.
5. If skin sensitisation or allergy occurs, stop using the product.
6. Do not store or expose to high temperatures or intense light.
7. To be used for dental purposes only.
8. Make sure the product is not dropped, especially onto the mandrel when using blocks.

WARRANTY

ITENA CLINICAL agrees to replace any product that is proven to be defective. ITENA CLINICAL does not accept responsibility for any damage or loss, direct or consequential, arising from the use of or inability to use the product. Before use, it is the user's responsibility to determine the suitability of the product for its intended use and the appropriate milling program for hybrid ceramic. The user assumes all risks and liability in that sense.

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