


# Dentoclic Scan Post

Geometry indicators required  
for scanning post-and-core cases



Enhance your  
daily practice

Compatible with Dentoclic

Simplified workflow

Precise impression



# Precise impression with a digital capture

## Key points

Compatible with Denclocic stainless steel posts and drills

Denclocic Scan Post library integrated into 3Shape™ solutions

Simplified workflow vs. traditional impression

Precise impression<sup>1</sup>: digital capture




Autoclavable



## Indications

- Dental prosthetic restorations
- Endodontic treatments

## Denclocic matching

| Scan Post  |       |             | Posts  |            | Drills  |                     |
|---|-------|-------------|---|------------|--|---------------------|
| Size  |       | Reference   | Size  | Reference  | Size   | Reference           |
| PL 1.2 (L 16 mm ø 1.2 mm)   | Long  | DESP-CCL1.2 | ○ L 11.5 mm ø 1.2 mm  | DIAL20-114 | ○ ø 1.2 mm   | DFA4-012 / FPPA-1.2 |
| PS 1.2 (L 12 mm ø 1.2 mm)   | Court | DESP-CCS1.2 | ○ L 9.6 mm ø 1.2 mm   | DIA20-095  | ○ ø 1.2 mm   | DFA4-012 / FPPA-1.2 |
| PL 1.3 (L 16 mm ø 1.3 mm)   | Long  | DESP-CCL1.3 | ● L 11.5 mm ø 1.3 mm  | DIJL20-114 | ● ø 1.3 mm   | DFJ4-095 / FPFJ-1.3 |
| PS 1.3 (L 12 mm ø 1.3 mm)   | Court | DESP-CCS1.3 | ● L 9.6 mm ø 1.3 mm   | DIJ20-095  | ● ø 1.3 mm   | DFJ4-095 / FPFJ-1.3 |
| PL 1.4 (L 16 mm ø 1.4 mm)   | Long  | DESP-CCL1.4 | ● L 13.5 mm ø 1.4 mm  | DIRL20-134 | ● ø 1.4 mm   | DFR4-115 / FPFJ-1.4 |
| PS 1.4 (L 12 mm ø 1.4 mm)   | Court | DESP-CCS1.4 | ● L 11.6 mm ø 1.4 mm  | DIR20-115  | ● ø 1.4 mm   | DFR4-115 / FPFJ-1.4 |
| PL 1.5 (L 16 mm ø 1.5 mm)   | Long  | DESP-CCL1.5 | ● L 13.6 mm ø 1.5 mm  | DIB20-135  | ● ø 1.5 mm   | DFB4-135 / FPFJ-1.5 |
| PS 1.5 (L 12 mm ø 1.5 mm)   | Court | DESP-CCS1.5 | ● L 13.6 mm ø 1.5 mm  | DIB20-135  | ● ø 1.5 mm   | DFB4-135 / FPFJ-1.5 |
| PL 1.6 (L 16 mm ø 1.6 mm)   | Long  | DESP-CCL1.6 | ● L 15.6 mm ø 1.6 mm  | DIV20-155  | ● ø 1.6 mm   | DFV4-155 / FPFV-1.6 |
| PS 1.6 (L 12 mm ø 1.6 mm)   | Court | DESP-CCS1.6 | ● L 15.6 mm ø 1.6 mm  | DIV20-155  | ● ø 1.6 mm   | DFV4-155 / FPFV-1.6 |
| PL 1.7 (L 16 mm ø 1.7 mm)   | Long  | DESP-CCL1.7 | ● L 17.6 mm ø 1.7 mm  | DIN20-175  | ● ø 1.7 mm   | DFN4-175 / FPFN-1.7 |
| PS 1.7 (L 12 mm ø 1.7 mm)   | Court | DESP-CCS1.7 | ● L 17.6 mm ø 1.7 mm  | DIN20-175  | ● ø 1.7 mm   | DFN4-175 / FPFN-1.7 |



[itena-clinical.com](http://itena-clinical.com)

Central Parc Bat B, 97 Allée de la Louve, 93420 Villepinte - FRANCE - T. + 33 1 45 91 61 40.  
FP.DENTOCICSCANPOST.EN-A-04/24

Enhance your daily practice

Follow us!



Medical device class I. Manufacturer : CoreMedTech ApS, Denmark. Reserved for professional dental use, not reimbursed by health insurance organizations under the LPPR. Carefully read the instructions and labeling before use.