

HydroSpeed HD¹²

Addition-crosslinking silicone-based
high precision impression material
(VPS*)



Convenient to use

Material saving

Compatible with all
impression techniques

Enhance your
daily practice



Hydrospeed HD

Addition-crosslinking silicone-based high precision impression material (VPS*)

Indications

Impression material

Detailed impression of inlay, onlay, crown and bridge preparation.

Bite registration

- Intra-oral bite registration
- Bite registration for orthodontic needs




Key points

Convenient to use

- > Easy impression taking, with elastic memory >99%
- > Thixotropic material (stays in place, does not sink)

29 to 38 % material savings thanks to tips

Compatible with all impression techniques



EXAMPLE OF IMPRESSION WITH HYDROSPEED HD

Impression material	
A. ● Light Body Quick 4 Cartridges 50 ml + 12 yellow mixing tips + 4 yellow intra-oral tips	HYLBQ4-50
B. ● Putty Soft Quick 1 catalyst pot 300 ml + 1 base pot 300 ml + 2 spoons	HYPQ2-300
Tips impression material	
C. 50 yellow mixing tips	MIXTYE-50
D. 50 yellow intraoral tips	IOTYE-50
Bite registration	
E. ● Bite Registration 2 cartridges 50 ml + 6 green mixing tips	HYREG2-50
Tips bite registration	
F. 50 green mixing tips	MIXTGR-50
G. 50 transparent flat tips	IOTTR-50
Accessorie ⁹	
50 ml 1:1 / 2:1 mixing gun	153241



itena-clinical.com

Central Parc Bat B, 97 Allée de la Louve, 93420 Villepinte - FRANCE - T. + 33 1 45 91 61 40.
FRHYDROSPEEDHD.EN-B-04/2024

Enhance your daily practice

Follow us!



* Polyvinylsiloxane.

¹² Medical device Class IIa. Manufacturer: Seil Global Co., Ltd, South Korea. Authorized Representative: MedNet GmbH, Germany.
⁹ Medical device Class I. Manufacturer: Medmix Switzerland AG, Switzerland. Authorized Representative: Medmix GmbH, Germany.
Reserved for professional dental use, not reimbursed by health insurance organizations under the LPPR. Carefully read the instructions and labeling before use.