

ISOsafe

Dental dams



Enhance your
daily practice

Hypersensitivity reduction

Tear and tensile resistance

Color contrast



Isolate the tooth, protect the patient

Key points

Hypersensitivity reduction

- > Powder concentration less than 2 mg
- > Non-latex IsoSafe: made of polyisoprene

Tear and tensile resistance

- > High elongation over 700%
- > Uniform thickness:
 - 0.28 mm IsoSafe latex
 - 0.25 mm IsoSafe non-latex
- > Tensile resistance:
 - ≥ 24 MPa IsoSafe latex
 - ≥ 20 MPa IsoSafe non-latex

Color contrast

- > Tooth visibility within the dental arch



IsoSafe latex dental dam

Size 6" x 6" (15,2 x 15,2 cm)
Average thickness (0,28 mm)



IsoSafe non-latex dental dam

Size 6" x 6" (15,2 x 15,2 cm)
Average thickness (0,25 mm)

Indication

Pre-cut latex and non-latex dental dams used as isolation barriers during endodontic and dental procedures.

IsoSafe

IsoSafe latex dental dam Box of 36 sheets	ISDALA6X6BL-36
IsoSafe non-latex dental dam Box of 15 sheets	ISDANL6X6PU-15



Enhance your daily practice

itena-clinical.com

Central Parc Bat B, 97 Allée de la Louve, 93420 Villepinte - FRANCE - T. + 33 1 45 91 61 40.
FPISOSAFE.EN | INDA 01/24

Follow us!

