

POLISH HD

Composite polishers



Convenience of use

Time saving

Natural aesthetic finish

Enhance your
daily practice



Fast & complete polishing

Key points

Convenience of use

- > **Twist:** polishes all surfaces with a single instrument. Medium-grit polishers (pre-polishing) and fine-grit polishers (polishing)
- > **Brush:** polishes the occlusal surfaces, including pits and fissures. Medium-grit polishers (pre-polishing) and fine-grit polishers (polishing)

Time-saving

- > 2 steps (instead of 4)
- > Use without polishing paste

Natural looking result

Reduced risk of microbial contamination:


- > Single use instrument




Simplified protocol



Twist
Anterior restorations



Brush
Posterior restorations



Also available



Cup
Lingual and palatal surfaces of anterior restorations, posterior cusps, proximal lines and angles.



Flame
Fissures, occlusal reliefs and cervical area.



Lens
Vestibular surface of incisors, composite surfaces, incisal angles and interproximal areas

Polish HD	
Set of 12 Flame polishers	PLFLAM-12
Set of 12 Cup polishers	PLCUP-12
Set of 12 Lens polishers	PLLENS-12
Set of 12 Brush polishers (medium-grit)	PLBRUSH-12
Set of 12 Brush polishers (fine-grit)	PLBRUSHF-12
Set of 12 Twist polishers (medium-grit)	PLTWM-12
Set of 12 Twist polishers (fine-grit)	PLTWF-12
Set of 6 assorted polishers (medium-grit) (x2 , x2 , x2)	PLASSORT-6
Set of 12 assorted polishers Medium-grit (x3 , x3) + Fine-grit (x3 , x3)	PLTATO-12
Set of 20 assorted polishers (medium-grit) (x5 , x5 , x5 , x5)	PLASSORT-20



Enhance your daily practice

itena-clinical.com

Central Parc Bat B, 97 Allée de la Louve, 93420 Villepinte - FRANCE - T. + 33 1 45 91 61 40.
FPPOLISHHD.EN | INDA 01/2024

Follow us!

