

# Provitemp

Short-term temporary  
cement (up to 6 weeks)



Enhance your  
daily practice

Tight seal

Easy handling

Suitable for use on living teeth



# Tight seal & eugenol-free

## Key points

### Tight seal


- > Low solubility ( $7\mu\text{g}/\text{mm}^3$ )

### Easy handling

- > Self-curing in 2 stages: gel for easy removal of excess and final set
- > Consistent, homogeneous mix for controlled application and distribution, thanks to the automix syringe and mixing tips

### Suitable for use on living teeth

- > Eugenol-free composition to avoid irritation and interference with the curing of resin-based definitive sealants



## Indications

Temporary cementation of crowns, bridges, inlays, onlays and splints.

### Provitemp

Provitemp Automix  
1x 6 g syringe + 10 mixing tips PTEMP1-10

Provitemp Automix  
2x 6 g syringes + 20 mixing tips PTEMP2-VP

### Tips

Refill with 20 mixing tips DTEM-20

Refill with 25 mixing tips + 25 thin intraoral tips  
(Length: 100 mm  $\varnothing$  1.1 mm) DCE-50

Refill with 25 mixing tips + 25 extra-thin intraoral tips  
(Length: 135 mm  $\varnothing$  0.09 mm) DCEXXF-50



[itena-clinical.com](http://itena-clinical.com)

Central Parc Bat B, 97 Allée de la Louve, 93420 Villepinte - FRANCE - T. + 33 1 45 91 61 40.  
FPPROVITEMPEN | IND.A 01/24

Follow us!

