

Traxodent

Hemostatic gingival
retraction paste



Effective retraction

Convenience of application,
rinsing and disposal

Hemostasis

Enhance your
daily practice



Retraction and hemostasis

Key points

Rapid, effective retraction in 2 min

- > Compression to stop bleeding while absorbing gingival fluids

Convenience of application, rinsing and disposal

Hemostasis

- > Spreading tissue in contact with air via chemo-mechanical expansion

Maneuverable, adjustable tips



Airtight resealable aluminum pouch

The paste does not dry out, and can be stored

Indication

Specifically formulated to provide chemo-mechanical gingival retraction and hemostasis. It is intended for use prior to taking an impression, cementation, cavity preparation or whenever hemostasis and retraction is required.

Traxodent Packs

2x 0,7 g syringes + 6 needle tips	TRASTRT-2
7x 0,7 g syringes + 15 needle tips	TRABOI-7
25x 0,7 g syringes + 50 needle tips	TRAPACK-25

Traxodent - Tips refills

20 needle tips 16 g	TRAEMB-20
60 needle tips 16 g	TRAEMB-60



itena-clinical.com

Central Parc Bat B, 97 Allée de la Louve, 93420 Villepinte - FRANCE - T. + 33 1 45 91 61 40.
FPTRAXODENT.EN | IND.A 01/24

Enhance your daily practice

Follow us!

